



Extraordinary vision



comes to life.



Experience SOLAOne™

Now available in 1.67 High Index and 1.67 Transitions®

Advanced design meets consumer insight
SOLAOne is based on the most comprehensive visual task analysis ever undertaken involving over 10,000 wearers.

Extraordinary performance, every day
SOLAOne provides exceptional performance across the 15 common visual task categories that make up daily life.*

Wearer's Guide CD-ROM
SOLAOne helps you create an extraordinary vision experience for your patients with an interactive wearer's guide CD-ROM – free with every pair.

Available with



Teflon
CLEAR COAT LENSES

Call 1-888-FIT-SOLA
or visit www.sola.com/fitSOLAOne/
to experience true versatility for yourself.

SOLA1One™
PROGRESSIVE LENSES

*Multi-task analysis and wearer trial data on file.

Now available in 1.67 High Index and 1.67 Transitions®

SOLA1One™
PROGRESSIVE LENSES




SOLAOne: Designed for real life.

Proven patient satisfaction

SOLAOne's success has been demonstrated in extensive clinical trials.

- Extremely low non-adapt rate, even when switching from other progressives.
- High 96% satisfaction rate among wearers.
- Consistently high satisfaction for visual tasks using every area of the lens.

SOLAOne Wearer Satisfaction Results*



*430 wearers recruited by the Queensland University of Technology School of Optometry and independent practitioners. Data on file.

Engineered for versatility

SOLAOne employs a unique, patent pending balance of central and dynamic vision for outstanding all-purpose performance.

- Generous distance with reduced peripheral blur and image swim
- Larger, wider intermediate vision area
- High, wide near vision area
- 18mm fitting height
- Advanced Design by Prescription™ technology

Available with




Teflon
CLEAR COAT LENSES

Available in a full range of premium lens options

SOLAOne HD		Rx Range	Add Powers
SOLAOne™ HD 1.67 High Index	Available Rx Teflon® coated only	-12.00 to +8.00 D	1.00 to 3.50
SOLAOne™ HD 1.67 Transitions® Gray	Available Rx Teflon® coated only	-12.00 to +8.00 D	1.00 to 3.50

SOLAOne Material Availability		Base Curves	Rx Range	Add Powers
NEW	1.67 MR-10 High Index	1.25, 2.75, 4.00, 5.75, 6.50, 7.50	-12.00 to +8.00 D	1.00 to 3.50
NEW	1.67 MR-10 High Index Transitions® Gray	1.25, 2.75, 4.00, 5.75, 6.50, 7.50	-12.00 to +8.00 D	1.00 to 3.50
	1.6 High Index	1.25, 3.25, 4.50, 6.25, 7.00, 8.00	-9.00 to +8.00 D	1.00 to 3.50
	Polycarbonate Transitions® Gray	1.25, 3.25, 4.50, 6.25, 7.00	-9.00 to +6.00 D	1.00 to 3.00
	Polycarbonate	1.25, 3.25, 4.50, 6.25, 7.00	-9.00 to +6.00 D	1.00 to 3.00
	Hard Resin Transitions® Gray	1.25, 3.25, 4.50, 6.25, 7.00	-9.00 to +6.00 D	1.00 to 3.50
	Hard Resin Transitions® Brown	1.25, 3.25, 4.50, 6.25, 7.00	-9.00 to +6.00 D	1.00 to 3.50
	Hard Resin	1.25, 3.25, 4.50, 6.25, 7.00	-9.00 to +6.00 D	1.00 to 3.50

Available from:




SOLA
www.sola.com